

Optum

Enterprise Imaging Amend a Report

March 24, 2025

This flyer will demonstrate how to amend a report within Ascend.

Amend a Report

From the patient's record:

STEP 1: Double-click the appropriate study to

open.

STEP 2: Click the **Ascend Report** icon.

1995

STEP 3: Click **Re-open**.

STEP 4: In the Re-open Amendment window, click

the Reason for Amendment text box.

NOTE: The amendment reason must be

provided to add an addendum text or edit report. This creates an accurate audit log and will not display on the

amended report.

Reason for amendment
Reason for amendment
Reason for amendment
Reason for amendment
Reason for amendment is information is not shown on the amended
report
Reason for amendment is required
Close
Close this window and leave the study for physician to amend

Addendum -- This information is shown on the amended report

Addendum -- This information is shown on the amended report

Sign Sign study without editing findings

Edit report Open the study for editing

Edit report Open the study for editing

NOTE: Without opening the report addendum text can be added in the Addendum box and use the

sign button to re-sign the report. This information is shown on the amended report.

STEP 5: Type the reason; then click **Edit report** to open and edit the report.

STEP 6: Add edits to the report.

STEP 7: When finished editing the report, click **Sign**.

STEP 8: Click **Confirm** in the Report Signature Confirmation window.



