

Optum

Enterprise Imaging Complete a Technologist's Cardiac Cath Ascend Report

March 5, 2025

This flyer will demonstrate how to begin, review, and assign a Cardiac Cath Ascend report.

Begin and Review Cardiac Cath Report

From the task list:

STEP 1: Click the Cardiac Cath study that needs to be reviewed and signed.


○○ Start reading ▼

- **STEP 2:** Click **Start reading** from the tool bar.
- **STEP 3**: From the **Startup** tab, click the type of study that was performed.
- STEP 4: Click Next.
- **STEP 5**: Select the appropriate **Case Type**.
- **STEP 6**: Click **Begin reporting**.
- **STEP 7**: Review the report and verify the information that has been

imported from the hemodynamic system is accurate.

STEP 8: Click **To be read** in the top toolbar.







Assign the Report

<u>STEP 1</u>: From the task list, right-click the **Sign icon** being sent to the cardiologist.



that is next to the task



- STEP 2: Click Assign.
- STEP 3: Click the Name field and start typing the cardiologist's name that the study will be assign to.
- **STEP 4**: Click the name being selected.
- <u>STEP 5</u>: Click the **Right Arrow** icon to move the selected name to the **Assign To** window.
- STEP 6: Click Assign.
- NOTE: Once the study has been assigned to the cardiologist,

the task will drop from the technologist's task list.

