

This flyer demonstrates how to Log in and out Enterprise Imaging Desktop.


Log into Enterprise Imaging (EI)

➤ From the desktop

STEP 1: Click the **Windows** icon at the bottom of the computer screen.



STEP 2: Click **Enterprise Imaging**.



STEP 3: Click the appropriate icon for logon to EI.

➤ Select appropriate desktop

- **Diagnostic** desktop is used by Cardiologists, Radiologists, Cardiology staff, Cath Lab staff, NUCs staff, and EMMC Vascular staff.
- **Acquisition** desktop is used by Technologists, Sonographers, and some PSRs.
- **Clinician** desktop is used by Providers, other Clinical roles, and office staff.



STEP 5: Type network **Username** and **Password**.




STEP 6: Click **OK**.

NOTE: The "Log on to" field should not be changed and should state ADMIN.

Log Off Enterprise Imaging

STEP 1: Click **File**.

STEP 2: Click **Log off**.



NOTE: If logging off the system, the EI login window will display to quickly log in when returning.

Exit Enterprise Imaging

STEP 1: Click **File**.

STEP 2: Click **Exit**.



NOTE: The system can be exited using the red X Close button in the top right corner. If the system is exited, access the system again through the Windows menu.