

This flyer demonstrates how to Log in and out Enterprise Imaging Desktop.

Log into Enterprise Imaging (EI)

➤ From the desktop


STEP 1: Click the **Windows** icon in the bottom left corner of the computer screen.



STEP 2: Click **Agfa Healthcare**.

STEP 3: Click **Enterprise Imaging**.

STEP 4: Click the appropriate icon for logon to EI.



➤ Select appropriate desktop

- **Diagnostic** desktop is used by Cardiologists, Radiologists, Cardiology staff, Cath Lab staff, NUCs staff, and EMMC Vascular staff.
- **Acquisition** desktop is used by Technologists, Sonographers, and some PSRs.
- **Clinician** desktop is used by Providers, other Clinical roles, and office staff.



STEP 5: Type network **Username** and **Password**.



STEP 6: Click **OK**.

NOTE: The "Log on to" field should not be changed and should state **ADMIN**.

Log Off Enterprise Imaging

STEP 1: Click **File**.

STEP 2: Click **Log off**.

NOTE: If logging off the system, the EI login window will display to quickly log in when returning.



Exit Enterprise Imaging

STEP 1: Click **File**.

STEP 2: Click **Exit**.

NOTE: The system can be exited using the red X Close button in the top right corner. If the system is exited, access the system again through the Windows menu.

