

From the Office of Health Informatics Enterprise Imaging (EI) Document and Image Management July 3, 2025

This flyer outlines how to Print, Import, Export, Fax, Email, Save to Media, and Transfer to DICOM exam documents from Enterprise Imaging (EI) Desktops.

Print, Fax, Save to Media Gallery, and Transfer DICOM data.

- **<u>STEP 1</u>**: Select **Patient**.
- **<u>STEP 2</u>**: Select **More** or **right-click**.
- **STEP 3:** Select **Print**.
- **<u>STEP 4</u>**: Select **Print and Export**.
- STEP 5: Select Distribute report, Save to media, or DICOM transfer.
- Distribute report -Print, Fax.
 - Print or Fax
 - Add local print
 - Select appropriate option under **Report sent to**.
 - Secure print, PDF option for xMedius faxing, Multi-Tech FaxFinder.
 - Select Distribute or Distribute & close.

				Print and expo	π		
	Import imag	es 🕨		Create print an	nd export task		
earches	Q QC task		9	Print label and	document		
criteria	PDF/Image	attachments	9	Print list			
T IIII	dy Add to My se	election	9	Save list			
t IDs	by Add to list				Report availab	le	
wailabla	Labels				Report availab	hidh	
a name	Download n	WO	-		Accession num	ohar	
e name				_	Accession num	ider	
5, 6 STAT! - I	Last updated at	09:44 č 2		All patient IDe	Accession num	Patient name and Patien	Mode
50	120/2022	uy date 🔺	B DE	S24544 - EMMC	NC202200004	O TESTING PROD2	NM
Report	NG, PROD1	ia DICOM tran	sfer	DICOM printing	1		
Report TESTIN	NG, PROD1 NG, PROD1	Distribute vi	a	DICOM printing	Report sent to	CC re,	
Report TESTIN	NG, PROD1 NG, PROD1 Name al print	Distribute vi	a a	DICOM printing	Report sent to	CC re • 1	8
Report TESTIN	IG, PROD1 Name al print	Distribute vi	a Ə	DICOM printing	Report sent to	d local print	8
Report TESTIN Cocc * default for Details for the	Name al print distribution the email	Distribute vi	a a	DICOM printing	Report sent to	d local print Add addressee	8
Report TESTIN Contemporta Contemport Contemport Contemport Contemport C	Name al print distribution the email	Distribute vi	a] 🖨	DICOM printing	Report sent to	d local print Add addressee	⊗
Report TESTIN Cocc * default for Details for t To Subject	Name al print distribution the email	Distribute vi	a 2015	DICOM printing	Report sent to	d local print Add addressee	8
Report TESTIN Coccentration * default for Details for t To Subject	Name al print distribution the email	Distribute vi	a 2015	DICOM printing	Report sent to	d local print Add addressee	8
Report TESTIN Cocc * default for Details for t To Subject Contents	IG, PROD1 IG, PROD1 IName If print Idistribution the email A report of	Distribute vi	a 2015	DICOM printing	r Report sent to INT	d local print Add addressee	8
5	IDs railable is name is . 6 STAT! - 	Add to My sr Add to list Add to l	Add to My selection Add to Ist Add to list Add to list Add to list Labels aliable Download now name s s s s Scheduled Study date Print and export	Add to My selection	Add to list Add to list Add to list Add to list Labels Download now aname s Add to list Add to lis	Add to My selection State ist Add to list Add to list Add to list Add to list Labels Download now Accession num Accession num Accession num Accession num Scheduled Study date All patient IDs Scheduled Study date All patient IDs Accession num. 09/30/2022 0624544 - EMMC NC20220004 Print and export	Add to My selection Study at 09:44 C aliable Download now Patient date of birth aname Accession number Assigned to (user+task assignment groups) 6 STATI - Last updated at 09:44 C Scheduled Study date All patient IDs Accession num Patient name and Patien 09/30/2022 Of Content of Con

Print and exp	ort	×
Distribute report	Save to media DICOM transfer DICOM printing	
🚹 Unable to d	istribute any reports; either there are no reports available or another user is still working on them.	
Report		

<u>NOTE</u>: "Unable to distribute any reports: either there are no reports available, or another user is still working on them," populates if nothing available to send at that time.

From the Office of Health Informatics EI – Document and Image Management July 3, 2025 Page 2 of 4

- Save to Media
 - Burn to CD/DVD
 - Select Driver.
 - Select Save to media or Save to media & close.
 - Exporting studies shows in the top banner beside Patient Name and MRN box.

Print and export						
Distribute report	Save to media	DICOM transfer	DICOM printing			
Select studies t	o export: 1 sele	cted				
EC Echo TESTINO	ocardiogram 2D G, CARD E (08/2	- 10/28/2022, 08:0 2/1984)	2 - 4775 images			
€ vậy						
Foi.nat	JPLG C BMF		OM h.e			~~~
Burn CD/DV Driver:	D -					
🗌 Include repo	rt as pdf					
				Save to media	Save to media & close	Close

<u>NOTE</u>: Once complete, a pop-up window appears stating: The following studies have been burnt to CD/DVD successfully with Name of study. Disk drive will open to retrieve the CD/DVD.

- > DICOM Transfer
 - Study is pre-selected.
 - Select **Destination** location from the dropdown.
 - Select **DICOM transfer** or **DICOM transfer & close**.
 - Exporting studies shows in the top banner beside Patient Name and MRN box.

Print and export	×
Distribute report Save to media V DICOM transfer DICOM printing	
Select studies to export: 1 selected	
EC Echocardiogram 2D - 10/28/2022, 08:02 - 4775 images TESTING, CARD E (08/22/1984)	
Study	
O Tagged images or series Reset	
Destination	-
NLH COMPASS - NLHCOMPASS	*
	_
DICOM transfer DICOM transfer & close Close	е

<u>NOTE</u>: Once complete, a pop-up window appears stating: The following studies have been transferred in accordance with filter configuration for the external system with Name of study.

<u>NOTE</u>: When exporting studies of the SAME Patient to CD or another destination, more than one can be bundled together.

From the Office of Health Informatics EI – Document and Image Management July 3, 2025 Page 3 of 4

Import CD

To import images to a selected study or a new order:

- **<u>STEP 1</u>**: Select the local drive under **In** the search area.
- **<u>STEP 2</u>**: Select study or studies to import.

<u>STEP 3</u>: Select **Import**.

🗢 Acquisition Desktop - 8.1.2 SP5 - ei81cs.a	agdataining net - El81 — _	
File Tools Listarea Help		
📙 Details 📰 Images 🕶 🌼 Open lis	ist 😰 Nee 💆 Import 🕶 Apre 🖛	
Started tasks (0)	Search 3	
	Search otheria DriD RW Drive (E) HIIH-202164838	
	Patient ID Patient's Name Patient's Birth Date	Y
	Study date (period) Study date Modatities in Study	٣
	1 study 🗛 📀	
	Patient ID Patient's Name Patient's Birth Study date Modalities In St Study description Accession num	
	970474244 O ⁷ 08-11-2018 CT ABDOMEN HMH202164938	

- **<u>STEP 4</u>**: Search for patient.
- **<u>STEP 5</u>**: Right-click in the patient row that needs a **QC task** to attach images to order entered.
- **<u>STEP 6</u>**: Hover over **QC task**.
- **<u>STEP 7</u>**: Left-click **Create QC task**.
- **<u>STEP 8</u>**: From the **QC action** dropdown, select the appropriate QC action, **Fix/Merge**.
- **<u>STEP 9</u>**: Enter comments, if necessary, in the **Comment** field.
- **<u>STEP 10</u>**: Click **Fix/Merge**.

Print Images

- **<u>STEP 1</u>**: Right-click an image.
- **<u>STEP 2</u>**: Select **Print and Export**.
- **<u>STEP 3</u>**: Select **Save to Media**.
- **<u>STEP 4</u>**: Select **Tagged images or series**.
- **<u>STEP 5</u>**: Select **Save to folder**.
- **<u>STEP 6</u>**: Select **JPEG**.

istribute report	Save to media DICOM transfer DICOM printing	
Select studies	to export: 1 selected	
CT Ab TESTI	Iomen and Pelvis with Contrast - 06/23/2014, 14:20 - 4 IG, PROD 13 L (10/10/1955)	4 images tagged
O Study		
O Key image	s Best v	
Tagged in	ages or series Reset	
Cours to fail	tor	
Folder	When emb or o/user data/EMMC/nswaf1/My Dou	
Format	JPEG BMP PNG DICOM file	
O Burn CD/D	VD	
Driver:	*	
Include rep	oort as pdf	

From the Office of Health Informatics EI – Document and Image Management July 3, 2025 Page 4 of 4

- **<u>STEP 7</u>**: Click **file folder** to browse location to save images.
- **<u>STEP 8</u>**: Select **My Documents**.
- STEP 9: Click Open.
- **<u>STEP 10</u>**: Go back to images, left-click the single image or series to tag.



- NOTE: Mindicates an image has been tagged.
- **<u>STEP 11</u>**: Go back to **Print and Export** screen after selecting all images to be printed.
- **<u>STEP 12</u>**: Select **Save to Media & Close**.

<u>NOTE</u>: Once complete, a pop-up window appears stating: The following studies have been exported to the local folder successfully with Name of study.

- **<u>STEP 13</u>**: Navigate to **This PC**, **Documents**.
- **<u>STEP 14</u>**: Double-click the **exported file**.
 - New window opens.
 - Double click the file.
- **<u>STEP 15</u>**: Select all images by clicking the first one, hold **Ctrl** key, and select the rest.
- **<u>STEP 16</u>**: Right-click one of the images.
- **STEP 17:** Select **Print**.
 - New Print Pictures window opens.
 - Select **Printer**.
 - Click **Print**.



For questions regarding process and/or policies, please contact your unit's Clinical Educator. For questions regarding workflow, please <u>place a ticket</u> to Health Informatics. For any other questions please contact the Customer Support Center at: 207-973-7728 or 1-888-827-7728.