
This Quick Reference Guide (QRG) demonstrates the customization of Enterprise Imaging.

Manage Activities Overviews

➤ From the EI desktop Collection – My list screen:

STEP 1: Click **List** area.

STEP 2: Click **Activities Overviews**.

- The **Manage Activities Overviews** window displays in the right pane.

STEP 3: Click the **Name** column header to sort the Overview list in alpha-numeric order.

NOTE: **Click-and-drag to expand/collapse columns, as desired.**

STEP 4: To set the default Overview, click the radio button of appropriate **Overview**.

NOTE: **The Clinician Workflow is the current default selection when first logging in to EI. There must always be a default Overview selected.**

STEP 5: Select the checkbox(es) of the Overview(s) in the **Show** dropdown list column to display in the Navigation area.

STEP 6: Click **Close**.

- The screen updates to display the selections.

Customize the Desktop

Customize the Task List

➤ From the EI desktop Task list – To do screen:

STEP 1: Click **List** area.

STEP 2: The **Customize columns**.

- The **Customize Columns** for the To Do window displays.

STEP 3: To change the order of the columns, click a column header in the **Displayed** field, then use the **Up** and **Down Arrows** to move it to the desired location.

NOTE: **The columns display left to right on the To Do list screen in the order they display in the Displayed field.**

STEP 4: To add a column header to view, click the column header in the **Available** field, then click the **Right Arrow** to move it the **Displayed** field.

NOTE: Any column header moved from the **Available** field to the **Displayed** field will automatically be at the bottom of the list; therefore, the last column displayed on the Task list. Use the **Up and Down Arrows** as shown to move a column to the desired location.

STEP 5: To remove a column header from view, click the column header in the **Displayed** field, then click the **Left Arrow** to move it to the **Available** field.

STEP 6: Click **OK**.

- The screen updates to reflect the changes.

Customize the Search Screen

➤ From the EI desktop Task list – To do screen:

STEP 1: Click the **Search** icon.

- The **Search** screen displays.

STEP 2: Click **List** area.

STEP 3: Click **Customize Search Fields**.

- The **Customize Search Fields** window displays.

STEP 4: Use the **Columns** options to select how many search field columns are in view.

STEP 5: To change the order of the search fields, click a search field in the **Displayed** search field list, then use the **Up** and **Down Arrows** to move it to the desired location.

NOTE: The search fields display left to right on the Search screen in the order they display is the **Displayed** search fields list and in the number of columns selected in the **Columns** option.

STEP 6: To add a search field, click the search field in the **Available** search fields list, then click the **Right Arrow** to move it to the **Displayed** search fields list.

NOTE: Any search field moved from the **Available** search fields list to the **Displayed** search fields list will automatically be at the bottom of the list; therefore, the last search field displayed. Use the **Up and Down Arrows** as shown to move a search field to the desired location.

STEP 7: To remove a search field from view, click the search field in the **Displayed** search fields list, then click the **Left Arrow** to move it to the **Available** search fields list.

- The search fields above update to display the changes immediately.

STEP 8: Click **OK**.

Use F12 Options

➤ From the EI desktop Task list – To do screen:

STEP 1: Click **Tools**.

STEP 2: Click **Options**.

NOTE: Press **F12** on the keyboard to open the Options window.

The Options window displays, open to the View tab of the Lists page.

These options have already been customized for everyone; however, they can be adjusted for user preferences. Use the “Restore default for Lists” button to automatically revert any options to the original default settings.

STEP 3: The **View** tab allows the lists and tasks to display expanded or collapsed, set the default display mode, choose the monitor set up preferences, open a study or task with a double-click, as well as other options.

STEP 4: The **Colors** tab represents STAT! and Urgent! tasks used throughout NLH.

STEP 5: The **Text** option allows customization of items related to the Reporting and Text area including view options, dictation controls, speech mike buttons, and foot pedals.

STEP 6: The **Images** option allows customization of items related to the Image area including image display options, keyboard shortcuts, and saving of mark-ups.

STEP 7: The **Other** option allows to turn on or off the Collaborator tool for the default login setting.

STEP 8: Click **Save** or **Save & Close** to save changes.