


This Quick Reference Guide (QRG) demonstrates how to complete a Physician Cardiac Cath report in Enterprise Imaging.

Complete a Physician Cardiac Cath Report


➤ From the Task List:

STEP 1: Click the study.


STEP 2: With your selected study highlighted, click **Start sign-off**.



STEP 3: Click **Start sign-off** again.



STEP 4: In the Ascend Reporting screen, click the Diagrams tab.



NOTE: When editing and making notes to the diagram on the left side, the changes will automatically appear in the Report on the right unless the Include in report check box is cleared. Include in report

Annotate the Diagram

STEP 1: Document the Coronary Dominance by clicking **Right**, **Left**, or **Co-**.



STEP 2: Click **Add/remove vessels**; then click the appropriate area to bring a vessel into view or remove it from view by graying it out.




STEP 3: Click the **X Close** button.

STEP 4: To document a lesion site, first select the lesion type; then click the area of the diagram where the lesion should display.

STEP 5: Document notes on the lesion in the pop-up menu; then, click the **X Close** button.

NOTE: To select a grayed out default item, click the Record Default Value icon. 

STEP 6: To document the lesion length on the Report side, click **Lesion length ?**. 

STEP 7: Select the appropriate number from the pop-up table or type the number in the field.

NOTE: Any field displaying a Question Mark denotes a suggested guideline by the American College of Cardiology (ACC).

To select a grayed out default item, click the Record Default Value icon. 

STEP 8: Click the **X Close** button.

STEP 9: To document that the patient had a prior stent, click **Prior stent**; then click the area of the diagram to place it.



STEP 10: Document more details in the pop-up Prior Interventions window. Then, click the **X Close** button.

NOTE: To select a grayed out default item, click the Record Default Value icon. 

STEP 11: Click the **Note Page** icon to add free-text notes to any line of the report. 

STEP 12: Click the text box and type note.

STEP 13: Click **Done**. 

NOTE: To delete a note, right-click and then select **Delete Finding** in the pop-up menu.

STEP 14: Once reviewed, the patient's findings are edited, and documentation is complete, click **Sign**.



STEP 15: Click **Confirm**. 