
This Quick Reference Guide (QRG) demonstrates how to amend a report in Enterprise Imaging.

Amend a Report

➤ **From the patient's record:**

STEP 1: Double-click the appropriate study to open.

STEP 2: Click the **Ascend Report** icon. 


STEP 3: Click **Re-open**. 

STEP 4: In the Re-open for Amendment window, click the **Reason for Amendment** text box.


NOTE: The amendment reason must always be documented to create an accurate audit trail; however, the reason entered here will not display on the amended report.


STEP 5: Type the reason; then click **Edit report** to open and edit the report. 

NOTE: Alternatively, you can click the **Addendum** text box to add documentation without fully opening the report.

STEP 6: Click **New summary item** to add new documentation to the Summary section. 

STEP 7: In the Add Note window, click the text box.

STEP 8: Type the note; then click **Done**. 

STEP 9: When finished editing the report, click **Sign**. 

STEP 10: Click **Confirm** in the Report Signature Confirmation window. 