
This Quick Reference Guide (QRG) demonstrates how to complete the Card Spect 4DM Phys workflow.

Complete the Card SPECT 4DM Phys Workflow

➤ From the Text area screen of a study:

STEP 1: Click the **Image area** icon.

- The Image area screen displays.



NOTE: The system may automatically launch in 4DM; if so, skip to Step 6. If 4DM does not launch automatically, navigate there using the Hanging Protocols icon.

STEP 2: Click the **Hanging protocols** icon.

- The Hanging protocols window displays.



STEP 3: Click the **Down Arrow** icon in the second field.

STEP 4: Click NM (Nuclear Medicine).

STEP 5: Click 4DM.

- The system displays in 4DM.

NOTE: The 4DM system is designed to perform its own manipulations and is often very accurate, requiring very little manual work to fix the images.

STEP 6: Click the **Polar Maps** tab.

- The area for calculating severity of defects displays.




STEP 7: Perform any edits or adjustments, as needed.

STEP 8: Click the **MPI Summary** tab.

- The screen provides a live display of the series slices with all the polar maps severity defects, as well as the calculated EF Curve.



NOTE: If this is a PET CT, click the **Reserve** tab, then review or adjust that section's post-processing.



STEP 9: Click the **Save** tab.

- The Save Data Status window displays.

STEP 10: Click OK.

- The saved documentation transfers to the Ascend report.
-

STEP 11: Click the **Clear Image** icon.

NOTE: At this point navigate to the Ascend report to sign it.

