

Lab values interfaced by Solarity require correction when incorrect.

## Incorrect Interfaced Lab Values

Uncharting incorrect lab values interfaced by Solarity can be done by Coding Operations.

### ➤ Locating the correct encounter

**STEP 1:** Search for practice specific encounter using the **Patient Search** tool.

- If already in the patient's chart, click the **LOC** link in the Banner Bar to select encounter.

**STEP 2:** Select the encounter.


**STEP 3:** Select **Yes** to navigate to the encounter.

- Assign relationship and select **OK** when prompted.

### ➤ Uncharting Interfaced Value


**STEP 1:** Navigate to dark side menu and click **Results Review**.

- Click **Labs** tab.




**STEP 2:** Change the timeframe if needed.

- Right-click in the gray bar and select **Change Search Criteria**.



**STEP 3:** Right-click the incorrect lab value and select **Unchart**. A prompt appears to enter a reason for uncharting.

- Uncharting a duplicate value, right-click the incorrect value and select **Unchart**.
  - Select the value to be uncharted and click **OK**.



**NOTE:** If multiple incorrect values, follow steps 1-3 to unchart each.

**NOTE:** Check to see if office staff entered correct values. If not, reach out to the end user to document the correct values using the Outside Labs PowerForm. Click [here](#) for flyer on Outside Lab Data Entry.