

On January 5, 2021, there will be changes to several Perioperative Doc segments within SurgiNet as a result of ongoing Periop Optimization Project work. These updates are necessary to steps toward upcoming upgrades. Removal of duplicative documentation, the addition of new functionality, and new information flow are among some of the exciting changes.

### What segments will be updated?

- Implant and Explant Log
- Skin Prep
- Hair Removal
- Counts Verification and Counts Action Taken
- Counts Verification – Additional
- Dressing/Packing, General Data, & Surgical Procedure Segments


### Implant and Explant Log Segments

The Explant Log segment will merge with Implant Log. Additionally, the Implant Log information will show up in the Implant History tool.

**STEP 1:** Within the IntraOp Record, select the **Implant Log** segment.


**STEP 2:** Select the procedure, then choose **Explant** when applicable.

- Explant Reason is on page 2.
- Click **Add** to Save.



**STEP 3:** Select the procedure, then choose **Implant**.

**STEP 4:** Begin typing and click the magnifying glass to search.



**STEP 5:** Select the implant from the list in the **Find Tool**, and click **Ok**.

- Manufacturer and Catalog numbers fill in for you, then move to page 2.

**STEP 6:** Lot/Batch number, Serial number, and Expiration date are entered here, as well as required Quantity.

- Click **Add** to Save, and both Explant and Implant are viewable for this procedure.

Procedure	Implant/Explant	Site	Site Comments	Description
Revision Arthroplasty ...	Explant			
Revision Arthroplasty ...	Implant	Knee		STRYKER KNEE INSERT ...

## Skin Prep

The Skin Prep segment is now consolidated into two pages, toggle between the two using the “1” and “2” found at the bottom of the segment list. Please note, the Hair Removal section will now be its own segment.

**STEP 1:** Within the Perioperative Doc tab, Select the **Skin Prep** section.

**STEP 2:** Select the appropriate prep agents, these are multiselect.

**STEP 3:** Select your prep site from the dropdown, select **>>** to add selection to the box.

**NOTE:** Do not select Add until both pages are complete for each Prep Site.

**STEP 4:** Continue to page 2.

**STEP 5:** Select the appropriate **Prep Area Detail** from the dropdown, select **>>** to add selection to the box.

**STEP 6:** Select **Prepped By** from dropdown. These names will populate based on those entered under **Case Attendees** segment. select **>>** to add selection to the box.

**STEP 7:** Complete remaining sections as appropriate.

**STEP 8:** Select **Add** when all sections are completed as appropriate for that Prep Site.

**STEP 9:** Repeat steps 2 through 8 for each additional Prep Site.

Prep Agents:

- 0.9% NaCl
- 5% Betadine Sol
- Alcohol 70%
- Baby Shampoo
- Betadine Scrub 7.5%
- Betadine Scrub 7.5% & Sol 10%
- Betadine Sol 10%
- Betadine Swab Stick
- BSS
- Chloraprep 2%
- Chlorhexidine Gluconate 4%
- Duraprep
- Exidine 4%
- Exidine 4% / Alcohol 70%
- Eye Lubricant
- Hibiclens
- Phisoderm
- Scrub Stat

Prep Site:

>> <<

Add Modify Remove Clear

Prep Area Details:

Prepped By:

Assjstive Device Used During Skin Preparation:

Immediate Skin Response to Prep:

○ No Skin Irritation  
○ Skin Irritation

Skin Prep Agent Dry Without Pooling:

○ Yes ○ No ○ N/A

## Hair Removal

The **Hair Removal** segment is a new addition to Perioperative Doc. This was part of the Skin Prep segment. This was also removed from “Safety Checklist 1 – Sign In” segment to reduce duplicative documentation.

**NOTE:** Sections with a Grey Box are required fields needing a selection.

**STEP 1:** Click one option in **Hair Removal Performed**.

**STEP 2:** Click one option in **Hair Removal Method**.

**STEP 3:** Click dropdown to select **Hair Removal Site**. Click to add one or more locations.

**STEP 4:** Click dropdown to select **Hair Removal Site Details**. Click to add one or more details.

**STEP 5:** Click dropdown to add **Hair Removed By**.

**STEP 6:** Click **Add**.

The screenshot displays the Perioperative Doc software interface. On the left, a navigation pane shows a list of documentation segments, with 'Hair Removal' highlighted in orange. The main window is titled 'Perioperative Doc' and shows a configuration screen for the 'Hair Removal' segment. The screen is divided into several sections: 'Performed' (with radio buttons for 'Preoperatively by Provider', 'At Home by Patient', 'In OR', and 'N/A'), 'Hair Removal Methods' (with checkboxes for 'Cliper', 'No hair removal performed', 'Patient performed own hair removal', 'Scissors', and 'Unable to determine method'), 'Hair Removal Site' (a dropdown menu), 'Hair Removal Site Details' (a dropdown menu), and 'Hair Removal By' (a dropdown menu). There are also 'Add', 'Modify', 'Remove', and 'Clear' buttons at the bottom. A 'Hair Removal Date/Time' field is visible on the right side.

## Counts Verification and Counts Action Taken

Counts Action Taken segment will be merged within the Counts Verification segment. Combining the two segments will ensure that required fields will not be missed.

- Added “**Procedure**” field to associate count performed to the appropriate procedure.
- Replaced Sponges with “**Soft Good**” and Other with “**Miscellaneous**”.
- Added a new field “**Interventions Taken Due to Missing Item**”.

The screenshot shows a form for 'Counts Verification and Counts Action Taken'. It has two main sections: 'Procedure' and 'Initial Counts'. The 'Procedure' section has a dropdown menu for 'Procedure:' with 'Arthroplasty Replacement Partial Hip' selected. The 'Initial Counts' section has a checkbox for 'Items included in the Initial Count:' and a list of items with checkboxes: 'Instruments', 'Sharps', 'Soft Goods', 'Miscellaneous', 'Throat Pack', and 'Bite block'.

- Updated Count Status to “Final Count Result” with updated options.

The image shows two parts of a software interface. On the left is a list of checkboxes under the heading "Interventions Taken Due to Missing Item:", including "Called for Assistance", "Intraoperative Radiograph of Surgical Site Obtained", "Room Searched", "Sterile Field Searched", "Supervisor Notified", and "Wound Reexamined". To the right is a dropdown menu titled "Final Count Result:" with the following options: "Correct (0.20 Met)", "Incorrect - Instruments (0.20 Not Met)", "Incorrect - Sharps (0.20 Not Met)", "Incorrect - Soft Goods (0.20 Not Met)", and "Incorrect - Miscellaneous (0.20 Not Met)". To the right of these is a section titled "Final Counts" containing a checkbox for "Surgical Team Notified of Count Results:" with radio button options for "Yes" and "No: See Segment Text".

- Added a new field “Added Surgical Team Notified of Count Results”.

### Counts Verification – Additional

This segment will be renamed to “Counts Verification – Cavity/Relief/Additional”. A new field will be added called “Count Time”.

The image shows a form field labeled "Count Time:". It consists of a text input box, a date/time mask "xx/xx/xxxx", and a dropdown arrow on the right side.

### Dressing/Packing, General Data, & Surgical Procedure Segments

The “Scope” question from General Case Data and “Closure Technique” from Dressing/Packing Segment will be moved to Surgical Procedures segment. It will be now named “Technique Details” group.

The image shows a section titled "Technique Details" with a checkbox. Below the checkbox are two radio button options: "Primary" and "Non-Primary". To the right of these is another checkbox labeled "Entire procedure was performed via laparoscope or robotic assistance:" with radio button options for "Yes" and "No".