


XERO Viewer is a medical image viewer that can be launched by an image link in the medical record or from a Web-Browser

Accessing XERO Viewer

- Open the patient's chart in PowerChart, navigate to the **Imaging component** on a **Workflow MPage**, and select the imaging study to be viewed.
 - Click the **View Image** icon. 

NOTE: Image link  in the EHR (Electronic Health Record) opens directly to the patient and study that is being opened.

Accessing XERO Viewer using a Browser


- Users not needing to use Ascend reporting.
 - Click to open link <https://brnagfaweb.me.emh.org/>.
- Cardiology and EMMC Vascular workflow to use Ascend reporting.
 - Click to open link <https://brnagfaweb.me.emh.org/?returnTo=/workflow/>.

NOTE: These links can be added to My Quick Links on the Intranet Home Page.

STEP 1: Enter **Patient Name**, **Date Of Birth**, **Accession number**, and **Study Date** (if known), and **Modality** to search for patient.


STEP 2: Click **Search**.

STEP 3: Select method of display, **Report and Images**, **Report Only**, or **Image Only**.




STEP 4: Double-click the study.

STEP 5: Select right arrow to show image area.



STEP 6: Select arrow on the bottom middle of the screen to scroll between images.




STEP 7: **Patient Banner** is used to verify patient.

STEP 8: Select the switch series icon in the upper right to switch between series.



STEP 9: Click one of the **Tools** to assist with reviewing images.



STEP 10: Select the ellipsis to change layout display.



STEP 11: Use the **Help** link to assist with any questions.

STEP 12: Click **Logout** once review of image(s) has been completed.

