

Using Related Results taskbar during Order Entry is an effective tool to review pertinent lab results at the time of entering orders.

Order Entry

Related Results via Orders profile can be collapsed and expanded depending on the need. The Related Results pane needs to be expanded to view results.

STEP 1: Enter Imaging contrast order from workflow MPage or Orders profile.

STEP 2: When Orders scratch pad opens, open **Related Results** pane to the left. For some positions, Related Results and Orders pane are collapsed. The expand arrow is available on the left.

STEP 3: Upon entering/signing order, review pertinent **Related Results**.

- For Imaging staff: **eGFR** and **Creatine Clearance** results details appear for all imaging contrast agents.

NOTE: Attachment of Order to Related Results is built individually by Contrast Agent. If new contrast is added to formulary, a request should be made to IS to have new contrasts added.

